

accessories

REFERENCE	ACCESSORIES	DESCRIPTION
7801	Stoppering Chamber	Glass cylindrical chamber with 3 non thermostated shelves and manual vial closing device. Includes 3 trays (for vials and serum)
8956	Chamber	Chamber 200 mm Ø, with 3 shelves and 3 trays (for vials, serum bottles and raw material)
7917	40 port manifold	40 tube manifold
70205	8 port manifold	Including 8 rubber three way valves CN 29/32 for flasks
7820	Adaptor plate for manifolds	Adaptor to place manifolds on the base unit



Stoppering chamber



Chamber



40 Port Manifold



8 Port Manifold

accessories technical data

	UNITS	CRYODOS -50	CRYODOS -80
Dimensions cylindrical chamber (ØxH)	mm	200 x 250	200 x 250
Tray dimensions (working Ø)	mm	160	160
Distance between shelves	mm	75	75

in-factory options

21OC2	Special vaccum pump	Vacuum pump specially prepared for freeze drying corrosive products
-------	---------------------	---

TELSTAR INDUSTRIAL, S.L.
SEDE CENTRAL
Josep Tapiolas, 120
08226 TERRASSA (España)
Tel. +34 937 36 16 00
Fax +34 937 85 93 42
E-mail: telstar@etelstar.com

TELSTAR USA, Inc.
Three Neshaminy Interplex
Suite, 301
TREVOSE, PA 19053 (USA)
Tel. +1 215 245 2280
Fax +1 215 245 2281
E-mail: broedel@telstarusa.com

DELEGACION MADRID
Santibáñez de Béjar, 2
28042 MADRID (España)
Tel. +34 913 71 75 25
Fax +34 913 71 79 41
E-mail: telstarm@etelstar.com

TELSTAR INT'L, Ltd.
Rm 712
Jingguang Centre Office Building
BEIJING (China)
P.O.: 100020
Tel. +86 10 65 26 11 13, 65261114
Fax +86 10 65 13 97 17
E-mail: telstar@telstar.com.cn



www.telstar.com